



Proposal

Visit us at <http://www.intech-ie.com>

Instrumentos y Equipos C.A

INTECH Instrumentos y Equipos C.A.

Av. Río Caura, Torre Humboldt, Ofic 04-04,
Caracas, Venezuela 1080
Tel: 58-212-9784430, 9784433, 977.7654, 977.6346
Fax: 58-212-977.3613
E-mail: sales@intech-ie.com
RIF J-30615033-1

Quote # C-2024-0084B
Date: 23/05/2024
Customer ID Sch - GUY

Bill to:
Guyana Logistics And Support Services Inc. (GLASS) Lot 87, Barrack Street, Kingston, Georgetown, Guyana. SA Phone: +592 2264090 Attn: Mohan Balram E-mail: mbalram@guyanalogistics.com

Ship to:
CIF Georgetown Airport. Custom Agent to be advised.

Sales Rep.	Your #	Our#	Destiny	Terms	Validity
N. Carricatti	Email	C-2024-0084/PL	Guyana	See below	30 days

Item	Cant	Description	Unit Price US\$	Total US\$
1	1	<p>ULTRAPYC 5000 PACK COMPLETE, ANTON PAAR PN PACK43704565</p> <p>This automatic gas pycnometer uses TruPyc technology for accurate results with sample amounts from 4.5 to 135 cm3 and its TruLock lid ensures repeatable results with a simple to use motion. These features and its Peltier temperature control make Ultrapyc 5000 instruments unmatched for true density measurement accuracy and repeatability. Add in its PowderProtect mode for safely measuring fine powders and Ultrapyc 5000 provides customers the widest range of gas pycnometer most complete gas pycnometer with the ability to accurately measure the widest range of sample sizes and types on a single instrument on the market today.</p> <p>Includes:</p> <p>*Standard UPYC, LARGE/MEDIUM/SMALL CELL WITH SLEEVE, CELL & SPHERES & Accessories Kit</p> <p>*REGULATOR (CGA)-N2/HE/AR/KR, 25 PSIG Two stage, gas regulator for inert gasses nitrogen, helium, argon, krypton, etc. North American CGA-580 fitting. Output pressure range 0 - 1.7 bar (0 - 25 psig). Includes isolation valve and 1/8 inch Swagelok output fitting.</p> <p>*UPYC, MICRO/MESO/NANO CELL WITH SLEEVE, CELL & SPHERES</p>	25,531.00	25,531.00
RECOMMENDED OPTIONALS				
2	10	<p>PYC, DISPOSABLE CUP (4 cm3), PACK OF 100 PN 243575</p> <p>These aluminum cups are intended to be used as disposable lightweight liners in the small pycnometer cell. Each cup weighs approximately 2 g and incorporates</p>	159.00	1,590.00

		a groove to facilitate easy removal using the standard pycnometer lift out tool. For use with viscous, sticky samples such as asphalt.		
		<p>Delivery time 4-5 weeks after receive your order and payment</p> <p>The delivery time of the items offered is subject to changes generated by events outside the company's control, such as:</p> <ol style="list-style-type: none"> 1. Climatic phenomena that do not allow the transport of raw material to the factory, local transport abroad of the manufactured material and/or international transport of the product. 2. Normal functioning of Venezuelan and/or foreign customs. 3. Disruption of public order in the country of origin, country of transit and/or final country. 4. Criminal and/or terrorist acts. 5. Any other event considered force majeure. <p>INCOTERMS CIF - Georgetown Airport</p> <p>Warranty <u>The proposal includes 3-year warranty. Warranty is valid as long as customer performs preventive maintenance after the first year of use. The warranty is extended for 1 year successively up to 3.</u></p> <p>Payment terms *100% Pre-payment with PO</p> <p>Volume & Weight Volume & Weight to be advised</p> <p>Travel Expenses Travel expenses for installation ARE NOT INCLUDED in this quotation. If special assistance is required for remote installation, it will be billed at 100.00 US\$ per working hour with a minimum of 8 hours.</p>		

Please Transfer to Account of:
Bank: Amerant Bank
Address: 220 Alhambra Circle, Coral Gables,
Florida 33134, USA
ABA: 067010509
Swift: MNBMUS33

Beneficiary: INTECH INSTRUMENTOS Y EQUIPOS
C.A.
Account # 8303790106

Subtotal	US\$	27,121.00
Freight	US\$	Included
Insurance	US\$	Included
Docs	US\$	--
TOTAL	US\$	27,121.00